



| | | | |
|---|---------------|----------|----------------|
|  | Record | N° | 520-FOR |
| | | Revision | 4 |
| | | Page | 1/12 |
| SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus | | | |

1. Introduction


The Summary of Safety and Performance (SSP) is one of the requirements of the new Regulation (IVDR 2017/746), specific for class C and D devices, to enhance transparency and adequate access to information. It intends to provide public access to summarised data on the safety and performance of class C and class D IVD devices to all intended users – professionals and lay persons.

2. Summary of Safety and Performance (SSP)

| Requirements based on IVDR Article 29 | Potential regulatory sources |
|--|---|
| Device identification and general information | |
| Name or trade name including any model number or version | FungaDia-Aspergillus Rapid Test |
| Manufacturer (name and address) | GaDia SA Route de l'Île-au-Bois 1A 1870 Monthey Switzerland |
| Manufacturers single registration number (SRN), if available | CH-MF-000031123 |
| Basic UDI-DI | 7649990065ASPMM |
| Intended purpose of the device | |
| Intended purpose and indications | FungaDia-Aspergillus is a rapid immunochromatographic test for the qualitative detection of the Aspergillus galactomannan antigen in serum and bronchoalveolar lavage (BAL) fluid from patients suspected of Fungal infections. This test is strictly for medical professional use only and not intended for personal use or home testing. The use of the test and the interpretation of the results should be done by a trained healthcare professional. The result of this test should not be the sole basis for the diagnosis; confirmation testing is required. |
| Target populations | Adult and pediatric patients with acute or subacute respiratory symptoms or fever or other suspicious symptoms or a known immunocompromised patient. The incidence of Invasive Aspergillosis (IA) in |

| | | | |
|---|---------------|----------|----------------|
|  | Record | N° | 520-FOR |
| | | Revision | 4 |
| | | Page | 2/12 |
| SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus | | | |

| | |
|---------------------------------|--|
| | immunosuppressed patients is rapidly increasing due to antibiotic abuse. |
| Contraindications (limitations) | <p>1. The product is only used for the detection of Aspergillus Galactomannan antigen in serum and BAL samples.</p> <p>2. The test results of this kit are for reference only and should not be used as the only basis for clinical diagnosis and treatment. The clinical management of patients should be comprehensively considered in conjunction with their symptoms, medical history, other laboratory tests and treatment responses.</p> |
| Device description | |
| Device description | <p>The principle of the test is colloidal gold immunochromatography. If the sample is positive, the antigens in the sample react with the red-colored nanoparticles and form a complex (Antigen - anti-Aspergillus monoclonal antibodies – gold nanoparticles), which was previously pre-dried on the conjugate pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-Aspergillus antibodies present on the membrane (Test line) capture the colored conjugate complex and a red line will appear. If the sample is negative, there is no Aspergillus antigens present or the antigens may be present in a concentration lower than the detection limit. The anti-Aspergillus antibodies present on the membrane (Test line) will not capture the antigen-red-colored conjugate complex (not formed), and the red line will not appear. Whether the sample is positive or not, the nanoparticle complex continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-mouse antibodies present on the membrane will react with the anti-Aspergillus antibodies coated on the gold nanoparticles and capture the complex to form a red line. The presence of this control red line serves as: (1) verification that sufficient volume is added, (2) that proper flow is obtained and</p> |

| | | | |
|---|---------------|----------|----------------|
|  | Record | N° | 520-FOR |
| | | Revision | 4 |
| | | Page | 3/12 |
| SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus | | | |

| | |
|--|---|
| | (3) an internal control for the reagents. The control line must always appear. |
| Reference to previous generation(s) or variants of the device (as applicable) and a description of the differences | N/A |
| Description of accessories intended to be used in combination with the device (as applicable) | N/A |
| Description of other devices and products intended to be used in combination with the device (as applicable) | Materials Required but Not Supplied 1. Pipettes and sterile tips 2. Timer 3. Disposable sterile micro-centrifuge tubes (eg. 72.692.005, Sarstedt) 4. Centrifuge 5. Heat incubator |

Standards Reference


| | |
|---|--|
| Harmonised standards and Common Specifications (CS) applied | IVDD 98/79/EC EN ISO 13485:2016 EN ISO 15223-1:2021 EN ISO 17511:2021 ISO 14971:2019 ISO 18113-1:2009 ISO 18113-2:2009 ISO 20417:2021 ISO 13975:2003 ISO 13612:2002 ISO 23640:2011 ISO 20916:2019 IEC 62366-1:2015+A1:2020 |
|---|--|

Summary of the Performance Evaluation

Methods:

A retrospective study was conducted at the University Hospital of Grenoble (CHU Grenoble-Alpes, Grenoble, France) using 153 serum samples and 33 samples from Bronchoalveolar lavage (BAL), retrospectively collected at the CHU Grenoble-Alpes in 2021 and 2022. 4 BAL samples (prevalence of 12%) and 49 serum samples (prevalence of 32%) were classified as positive with the CE-IVD PLATELIA ELISA Aspergillus Ag Galactomannan assay (BioRad, Marne-la-Coquette, France), according to manufacturer instruction for use.

To perform the test, 0.3 ml of defrosted serum sample was diluted in 0.1 ml of treatment buffer, incubated at 130°C for 6 min and centrifuged 10 min at 10'000xg. For BAL samples, 0.3 ml of samples were heated at 100°C for 3 min and centrifuged 10 min at 10'000xg. Then, 50 ul of

| | | | |
|---|---------------|----------|----------------|
|  | Record | N° | 520-FOR |
| | | Revision | 4 |
| | | Page | 4/12 |
| SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus | | | |

treated sample were added to the sample hole and incubated at room temperature for 10 minutes. The result was read visually and in parallel with the LFA Reader (Qiagen) to determine the signal intensity (in mV) of test line, blinded to the reference ELISA result. The signal intensity and results of the test were reported and analyzed.

The primary end point was to assess the diagnostic performance of FungaDia-Aspergillus Rapid test (GaDia SA, Monthey, Switzerland) in serum and BAL samples against the ELISA reference method PLATELIATM ELISA Aspergillus Ag Galactomannan assay (BioRad, Marne-la-Coquette, France)

Vassarstats online tool (www.vassarstats.net) was used to calculate sensitivity (SE), specificity (SP), positive and negative predictive values (PPV, NPV), 95% confidence intervals, median, and Interquartile range (IQR); while significance (p-values) was calculated using student t test for independent samples with equal variances. Statistical significance was defined as $p < 0.05$. ROC curve analysis was performed using easyROC (<http://www.biosoft.hacettepe.edu.tr/easyROC/>). Box plot are generated using BoxPlotR (<http://shiny.chemgrid.org/boxplotr/>).

Results:

The diagnostic performance of FungaDia-Aspergillus rapid test on serum samples is described below in the Table.

| Platelia (BioRad)/Clinical diagnostic | | | |
|---------------------------------------|-------|---------------------|-----|
| | | + | - |
| FungaDia ICA | + | 38 | 4 |
| | - | 11 | 100 |
| Sensitivity: | 77,6% | (CI95%: 63.0-87.8%) | |
| Specificity: | 96,2% | (CI95%: 89.9-98.8%) | |
| PPV: | 90,5% | (CI95%: 76.5-96.9%) | |
| NPV: | 90,1% | (CI95%: 82.6-94.7%) | |

Table 2: Diagnostic performance of FungaDia-Aspergillus in Serum samples compared to ELISA Platelia

A total of 11 discordant results showed a negative FungaDia-Aspergillus results while positive with ELISA reference method (false negative results) and 4 false-positive results.

The kappa coefficient was calculated between FungaDia-Aspergillus and reference ELISA method for Serum samples. A kappa coefficient of 0.77 was obtained

The diagnostic performance of FungaDia-Aspergillus rapid test on BAL samples is described below in the Table. Only the results of 30 BAL samples were collected. 1 sample was invalid (no flow, tracheal aspiration sample) and 2 were not tested due to low number of rapid test kit.

SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus

Platelia (BioRad)/Clinical diagnostic

| | | + | - |
|--------------|--------|---------------------|----|
| FungaDia ICA | + | 2 | 1 |
| | - | 0 | 27 |
| Sensitivity: | 100,0% | (CI95%: 20.0-100%) | |
| Specificity: | 96,4% | (CI95%: 79.8-99.8%) | |
| PPV: | 66,7% | (CI95%: 12.5-98.2%) | |
| NPV: | 100,0% | (CI95%: 84.5-100%) | |

The sensitivity of FungaDia-Aspergillus is good on BAL samples due to low number of samples available and the wrong sample preparation used for the half of samples. These samples were prepared according to BioRad sample preparation.

The specificity of the rapid test is good and similar to Serum samples.

Due to low number of positive BAL samples, the sensitivity of the test need to be interpreted with caution. The kappa coefficient is 0.78 for BAL samples.

Discussion & conclusion:

The main finding of this evaluation study, using an unmatched case control design including 32% (49/153) of positive samples for Serum samples and 12% (4/34) positive samples in BAL positive group, is that the diagnostic accuracy of FungaDia-Aspergillus RDT on serum samples when compared to ELISA Platelia confirmed cases displayed a SE of 78 %, a SP of 96%, a PPV of 91 % and a NPV of 90%, respectively. For BAL samples, due to the low number of positive cases, the performance values need to be evaluated with care. The SE is 100%, a SP of 96 %, a PPV of 67 % and a NPV of 100 % for BAL samples.


Aspergillus galactomannan detection in serum and BAL samples are widely used for the diagnosis of Invasive Pulmonary Aspergillosis (IPA). Most of the actual procedure are done in laboratory environment such as ELISA. D'Haese et al (2012) have assessed the clinical validity of such detection. Their conclusions are that detection of GM in BAL fluid samples of patients at risk of IPA has an excellent diagnostic accuracy. Recently, several guidelines have demonstrated the utility of GM detection in BAL fluid. Less data is available for serum samples and the recommendation are for specific cases only.

Gupta et al (2017) have showed that detection of Aspergillus galactomannan in BAL is more sensitive (87.5%) compared to the detection in serum samples (45%). Similarly, other recent studies confirm the observations (Park et al. 2011; Sehgal et al. 2019; Wu et al. 2021).

In order to compare the agreement between FungaDia-Aspergillus rapid test and ELISA, the kappa coefficient was used (Viera & Garrett 2005). When assessing the ability of a test to be helpful to clinicians, it is important that its interpretation is not due to chance. The kappa coefficient gives a valuable information regarding the agreement between 2 observers or tests.

FungaDia-Aspergillus has a kappa coefficient of 0.77 when comparing serum samples with ELISA method and 0.78 for BAL sample. This coefficient is defined as "Moderate agreement" (Viera & Garrett 2005).

There are however some limitations to this study. First, this is the result of a method evaluation study and not a seroprevalence study. Therefore, the PPV obtained here (based on a 25-30 %

| | | | |
|---|---------------|----------|----------------|
|  | Record | N° | 520-FOR |
| | | Revision | 4 |
| | | Page | 6/12 |
| SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus | | | |

prevalence) will be lower in a low prevalence setting. The test was performed in a laboratory environment.

In conclusion, FungaDia-Aspergillus is not meant to replace classical laboratory diagnostic such as ELISA or EIA but could be considered in low-income countries or where laboratory equipment is not available. It could be used as a pre-screening tool to detect Aspergillus galactomannan in Serum and start the appropriate treatment earlier

In conclusion, FungaDia-Aspergillus gives valuable information in only 15 minutes to start the treatment earlier. Performance with ELISA methods, longer and requiring laboratory equipment, has moderate agreement with Serum Samples. A larger study with Bal samples should be conducted.

Summary of the Post-Market Performance Follow-Up (PMPF)

Two PMPF were performed in 2022-2023, one at the University Hospital of Besançon (France) and the second one at the University Hospital of Grenoble (France).

Comparison of different lateral flow assays on bronchoalveolar lavage fluid for invasive aspergillosis screening in non-hematological patients

Devillers Eliane¹, Scherer Emeline^{1,2}, Navellou Jean-Christophe³, Grenouillet Frédéric^{2,4}, Millon Laurence^{1,2}, Bellanger Anne-Pauline^{1,2}

¹Parasitology-Mycology Department, University Hospital of Besançon, 25000 Besançon, France

²CNRS-University of Franche-Comte/ UMR 6249 Chrono-environment, 25000 Besançon, France

³Intensive Medical Care Department, University Hospital of Besançon, 25000 Besançon, France

⁴Infectious Serology Department, University Hospital of Besançon, 25000 Besançon, France

Several lateral flow assays (LFA) capable of detecting Aspergillus fumigatus in serum and broncho-alveolar lavage fluid (BALF) within the hour, thereby potentially accelerating the screening process, are now commercially available. We prospectively compared three LFA targeting A. fumigatus on BALF collected from non-surgical intensive care patients between June 2022 and February 2023 at the University Hospital of Besançon (France). The three LFA tested were Sōna Aspergillus galactomannan LFA (Immy), Fungadia Aspergillus antigen (Gadia), and AspLFD (OLM Diagnostics). We compared the results of these LFA with those of the galactomannan (GM) Platelia Aspergillus enzyme immunoassay (Bio-Rad), culture on Sabouraud medium and Aspergillus qPCR.

From June 2022 to February 2023, we tested BALF samples from ICU patients with suspected IMD with the three LFA, in accordance with the manufacturers' instructions. The results were read without a digital reader. The results of the three LFA tested were compared with those for Platelia GM detection, culture on Sabouraud medium for two weeks at 30°C and 35°C and the results of Aspergillus qPCR. GM detection with the Platelia Aspergillus immunoassay (Biorad, Marnes-la-Coquette, France) assay was performed with the cutoff suggested by the manufacturer (index of 0.5) used to classify samples as positive or negative. For Aspergillus qPCR, we used an in-house assay combining a mitochondrial and a ribosomal target, considered negative if the cycle of quantification (Cq) was >45.

We tested 97 BALF samples from 92 ICU patients with a median age of 61 years [18-88 years] prospectively. These 97 BALF samples included 35 hemorrhagic samples (36%).

Overall, 84 of the 97 BALF samples tested negative in all three LFA studied (IMMY, OLM and GADIA); these BALF samples also tested negative for Aspergillus species by culture on Sabouraud medium, GM detection in the Platelia assay and our in-house Aspergillus qPCR.

Four of the BALF samples tested positive only in the OLM assay, giving negative results in all the other tests (IMMY, GADIA, Platelia, culture and Aspergillus qPCR) (Table 1).

Only one BALF sample tested positive in all three LFA assays, as well as the Platelia assay (index of 3.5). However, this sample tested negative by in-house Aspergillus qPCR and by culture on Sabouraud medium, which was instead positive for *Candida albicans* (Table 1).


Three BALF samples tested negative in all three LFA, but positive for GM detection with the Platelia assay (index of 0.58, 0.66 and 0.81). For these samples, culture on Sabouraud medium and in-house Aspergillus qPCR were negative (Table 1).

Five BALF samples tested negative with the three LFA and in the Platelia assay (index of 0.06, 0.18, 0.26, 0.37 and 0.41) but positive by in-house Aspergillus qPCR (Cq 27,32,35,37 and 39). Two of these samples tested positive for *A. fumigatus* by culture on Sabouraud medium (Table 1).

Table 1: Summary of the results for the 97 BALF samples from 92 patients tested

| | LFA Immy | LFD OLM | FUNGADIA Gadia | Galactomannan Platelia | Culture on Sabouraud medium | Aspergillus qPCR |
|-----------------------|----------|----------|----------------|---|---|------------------------|
| 84 BALF (80 patients) | Negative | Negative | Negative | Negative | Negative | Negative |
| 4 BALF (4 patients) | Negative | Positive | Negative | Negative (GM index 0.03,0.04,0.08,0.22) | Negative | Negative |
| 1 BALF (1 patient) | Positive | Positive | Positive | Positive (GM index >3.5) | Positive > 1000 CFU/mL <i>C. albicans</i> | Negative |
| 3 BALF (3 patients) | Negative | Negative | Negative | Positive (GM index 0.58,0.66,0.81) | Negative | Negative |
| 3 BALF (2 patients) | Negative | Negative | Negative | Negative (0.06,0.37,0.41) | Negative | Positive (Cq 32,35,39) |
| 2 BALF (2 patients) | Negative | Negative | Negative | Negative (0.18,0.26) | Positive <i>A. fumigatus</i> | Positive (Cq 27,37) |

This study at CHU Besançon was presented on a poster during the French Congress of Medical Mycology.



MARRAKÉCH 24-27 Mai 2023

E. Devillers¹, E. Scherer¹, F. Grenouillet², J.C. Navellou³, L. Millon¹, A.P. Bellanger¹

1.Parasitologie-Mycologie, CHU de Besançon, Besançon, France; 2.Sérologie fongique et parasitaire CHU de Besançon, Besançon, France; 3.Réanimation médicale, CHU de Besançon, Besançon, France;

INTRODUCTION

Afin d'améliorer la rapidité du diagnostic des aspergilloses invasives, plusieurs tests immuno-chromatographiques permettant de détecter la présence d'*Aspergillus fumigatus* dans l'heure, sur sérum et sur lavage broncho-alvéolaires (LBA) ont été mis sur le marché ces dernières années (1-3).

MATERIELS ET METHODES

Dans cette étude, 3 coffrets ont été évalués prospectivement de juin 2022 à février 2023, pour la détection d'*A. fumigatus* dans des LBA de patients de réanimation médicale:

- ✓ Aspergillus GM LFA (IMMY)
- ✓ AspLFD (OLM)
- ✓ FungaDia Aspergillus antigen (GADIA)

Un LBA positif en galactomannane (GM) Platelia Biorad, en culture et en qPCR Aspergillus (Cq 36, technique maison) a servi à valider les 3 coffrets (Figure 1).

Les résultats obtenus ont été comparés à ceux obtenus en Platelia, à la culture sur milieu Sabouraud en 14 jours ainsi qu'à la qPCR Aspergillus (4).



Figure 1

RESULTATS

97 LBA de 92 patients différents ont pu être testés avec les 3 coffrets:

- ✓ 84 LBA étaient négatifs avec les 3 tests
- ✓ 4 LBA étaient positifs uniquement avec le kit AspLFD (OLM) (Tableau 1)
- ✓ 1 LBA était positif avec les 3 tests: le GM Platelia était aussi positif (>3,5) mais culture négative pour *A. fumigatus* (positive pour *Candida albicans*) et la qPCR Aspergillus était négative.

Par ailleurs sur cette série,

- ✓ 3 LBA étaient positifs uniquement en GM Platelia (Tableau 1)
- ✓ 5 LBA étaient positifs en qPCR Aspergillus dont 2 positifs également en culture pour *A. fumigatus* (Tableau 1)

Sur 92 patients testés,

- ✓ 1 diagnostic d'aspergillose possible
- ✓ 1 diagnostic d'aspergillose probable

| | LFA Immy | LFD OLM | FungaDia Gadia | GM Platelia | Culture | qPCR Aspergillus |
|--------|----------|---------|----------------|-------------------------------|-----------------------------|-----------------------|
| 84 LBA | négatif | négatif | négatif | négatif | négatif | négatif |
| 4 LBA | négatif | positif | négatif | Négatif (0.03,0.04,0.08,0.22) | négatif | négatif |
| 1 LBA | positif | positif | positif | Positif (>3.5) | Positif <i>C. albicans</i> | négatif |
| 3 LBA | négatif | négatif | négatif | Positif (0.58,0.66,0.81) | négatif | négatif |
| 3 LBA | négatif | négatif | négatif | Négatif (0.06,0.37,0.41) | négatif | Positif (Cq 32,35,39) |
| 2 LBA | négatif | négatif | négatif | Négatif (0.18,0.26) | Positif <i>A. fumigatus</i> | Positif (Cq 27,37) |

Tableau 1

DISCUSSION

Cet essai prospectif est le seul à inclure le test FungaDia Aspergillus antigen (GADIA). Les deux autres tests, Aspergillus GM LFA (IMMY) et AspLFD (OLM) ont déjà été évalués dans d'autres études (1-3).


Dans cet essai, nous avons malheureusement eu peu de prélèvements positifs.

Nous avons cependant pu évaluer la praticabilité de chaque coffret. Le kit AspLFD (OLM), dont le protocole est le plus simple, a présenté le plus de faux positifs (4/97); la bande est également moins nette en cas de positivité (Figure 1).

Le test FungaDia Aspergillus antigen (GADIA) s'est positif après le temps d'incubation recommandé par le fournisseur dans 6 cas.

Le test Aspergillus GM LFA (IMMY) présente quant à lui l'inconvénient d'un conditionnement par 50 bandelettes avec une date de péremption courte par rapport aux autres coffrets (20-25 tests / coffret et date de péremption plus longue).

Références: 1. Scharmann et al, 2020, Mycoses; 2. Mercier et al, 2020, Critical Care; 3. Aerts et al, 2022, J Clinical Mycology; 4. Bellanger et al, 2015, Med Mycol

| | | | |
|---|---------------|----------|----------------|
|  | Record | N° | 520-FOR |
| | | Revision | 4 |
| | | Page | 8/12 |
| SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus | | | |

Retrospective & prospective evaluation using LBA and serum from adult patients at CHU Grenoble-Alpes.

Muriel Cornet¹ & Danièle Maubon¹

¹ Parasitology-Mycolology Laboratory, Institut de Biologie et de Pathologie, CHU de Grenoble-Alpes (France).

Aims of the study are:

- 1.1. To evaluate the contribution of the rapid test for the detection in the context of emergency of galactomannan in LBA/serum:
 - in practice: we could at each series of antigen assay, check the positives with the RDT in order to establish a detection threshold (test sensitivity?).
Propose an insightful positioning of the test in the diagnostic tree
- 1.2. Work on the grey areas of the tests: on samples very slightly < 0.5 or very slightly > 0.5 (e.g. between 0.4 and 0.6?) True positives? false negatives? in this problematic we will also extract data on serum passed 'in kinetics' (with one or more priors)
- 1.3. Work on cross-reactions (especially interferences with immunoglobulins): in practice: a study of this type has already been done for BDG and platelia.
[https://www.clinicalmicrobiologyandinfection.com/article/S1198-743X\(20\)30140-3/fulltext#tbl1](https://www.clinicalmicrobiologyandinfection.com/article/S1198-743X(20)30140-3/fulltext#tbl1)

Methods:

A total of 563 samples were found positive with Ag Platelia > 0,5 in the biobank of CHU Grenoble-Alpes from 2017 to 2022.

According to Aspergillus team of CHU Grenoble, 306 samples were considered as True positive (TP), corresponding to 109 patients, and 257 samples were considered as False positive (FP), corresponding to 42 patients. The figure below described the tested population and exclusion criteria.


For this evaluation, 24 samples FP and 24 samples TP, retrospectively collected, were tested in parallel with rapid test kit FungaDia-Aspergillus.

Prospective samples collected on December 2022 (n= 38) were also evaluated with the rapid test blinded to the reference test Platelia ELISA kit (BioRad).

The retrospective clinical evaluation was conducted at the University Hospital of Grenoble in December 2022 on 46 characterized and archived serum samples (figure 1).

**SUMMARY OF SAFETY AND PERFORMANCE (SSP)
FungaDia-Aspergillus**

Figure 1:
Flow Chart



COVID19 Associated Pulmonary Aspergillosis

These samples were collected between 2017 and 2022 and included CAPA (COVID-Associated Pulmonary Aspergillosis) patients. Only True positive samples and False positive samples on Platelia were used for this study, in order to exclude any bias in the evaluation. The FungaDia-Aspergillus tests were performed according to the manufacturer's instructions. Prospective samples n=38 (9 BAL and 29 serum samples) were also tested on December 2022 using the rapid test. Sensitivity, specificity, positive and negative predictive values (PPV, NPV) using the Platelia ELISA kit (Bio-Rad) results as reference were calculated.

Results:


Retrospective study (n=46)

The overall performance of Rapid test kit FungaDia is presented in the figure below. A total of 46 positive samples were used for the retrospective evaluation of the rapid test. The overall sensitivity of the kit is 89%. As no negative samples were evaluated, the specificity was not evaluated in this retrospective study.

Retrospective study

| FungaDia RDT | Platelia | |
|--------------|----------|---|
| | + | - |
| + | 41 | 0 |
| - | 5 | 0 |

Sensitivity: 89,1% (CI95%: 75,6-95.9%)

| | | | |
|---|---------------|----------|----------------|
|  | Record | N° | 520-FOR |
| | | Revision | 4 |
| | | Page | 10/12 |
| SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus | | | |

More negative samples will be tested in the prospective study presented below. The performance of the kit based on the retrospective evaluation is better than the performance tested in January 2022 (sensitivity of 78%)

Prospective study (n=38)

The prospective evaluation was conducted in December 2022, directly after the evaluation with the reference test (Platelia, BioRad) blinded to the reference test. A total of 10 positive samples (6 BAL, 4 serum) and 28 negative samples (3 BAL and 25 serum) were tested. The performance of the kit is presented below.

The sensitivity of the kit was 70% and the specificity 89%. These results must be evaluated with care as only 10 positive samples were evaluated.

| Prospective Study | | |
|---------------------|--------------|---------------------|
| | Platelia | |
| FungaDia RDT | + | - |
| + | 7 | 3 |
| - | 3 | 25 |
| Sensitivity: | 70,0% | (CI95%: 35,4-91.9%) |
| Specificity: | 89,3% | (CI95%: 70,6-97.2%) |
| PPV: | 70,0% | (CI95%: 35,4-91.9%) |
| NPV: | 89,3% | (CI95%: 70,6-97.2%) |

Overall performance (n=84)


The overall performance of the kit, combining prospective and retrospective samples (n=84) is presented below on the table. The overall sensitivity is 85.7% and the specificity 89.3%.

| | Platelia | |
|---------------------|--------------|---------------------|
| FungaDia RDT | + | - |
| + | 48 | 3 |
| - | 8 | 25 |
| Sensitivity: | 85,7% | (CI95%: 73,2-93.2%) |
| Specificity: | 89,3% | (CI95%: 70,6-97.2%) |
| PPV: | 94,1% | (CI95%: 82,8-98.5%) |
| NPV: | 75,8% | (CI95%: 57,4-88.3%) |


Conclusion:

The performance of rapid test is comparable with similar product on the market with an overall performance (2 studies combined) on serum and BAL samples (n=267) of a sensitivity of 82% and a specificity of 95%.

| Metrological traceability | |
|--|-----|
| Metrological traceability of assigned values | N/A |

| | | | |
|---|---------------|----------|----------------|
|  | Record | N° | 520-FOR |
| | | Revision | 4 |
| | | Page | 11/12 |
| SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus | | | |

| Users | |
|--|--|
| User Profile | The tests can be performed in laboratories by health care workers or laboratory technicians with appropriate training in sample collection, biosafety and in the use of rapid tests. |
| User Training | Appropriate training in sample collection, biosafety and in the use of rapid tests. |
| Device Risks Information | |
| Residual risks and undesirable effects | <ul style="list-style-type: none"> - Contamination of the user by infected samples - Wrong interpretation of the test results - False negative - Interference - Cross-reactivity |
| Warnings and precautions | <ol style="list-style-type: none"> 1. This product is used for in vitro diagnosis, professional use only. 2. Do not reuse the test. Do not use the test after expiry date 3. Please read the test results within the specific time to avoid wrong interpretation. 4. Do not use the components from different batches or different types of reagents. 5. Properly dispose the specimen and used materials following the local biohazardous disposal regulation. 6. Use protective equipment when handling samples and tests as they may contain infectious agents and human or animal components. 7. Sodium azide is used as preservative in the sample treatment solution. Dispose material according to relevant local regulations and avoid contact with eyes and skin. 8. When the content of Aspergillus antigen in the sample is very high, the line C may be weakened. 9. Very high concentrations of Aspergillus antigen cause a hook-like effect, leading to false negative results. In this case, it is recommended to use a physiological saline |

| | | | |
|---|---------------|----------|----------------|
|  | Record | N° | 520-FOR |
| | | Revision | 4 |
| | | Page | 12/12 |
| SUMMARY OF SAFETY AND PERFORMANCE (SSP) FungaDia-Aspergillus | | | |

| | |
|--|--|
| | <p>solution to make a 5-10 times dilution of the sample.</p> <p>10. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established</p> |
|--|--|